



Dentists Recommend Bristle Tip Brushing

This is a brushing technique that involves brushing the teeth with the tip of the brush, using soft strokes. By maneuvering the bristles gently, plaque is removed more effectively, while being gentle to the gums.



Gently brush using the tips of the bristles.



Using so much force such that the bristles bend when brushing.



DR.SATO
The Nippon Dental University, Japan
School of Life Dentistry of Niigata

Bristle Tip Brushing is effective to maintain the health of our teeth.

Even if one were to apply much strength when brushing, it is difficult to get rid of the plaque. On the contrary, it will hurt the gums. We dentists recommend, brushing by using the tips of the bristles, also known as 'Bristle Tip Brushing'. Plaque can be removed more effectively by moving the brush on the teeth numerous times than by making large brushing movements. Maintain the health of your teeth and gums using 'Bristle Tip Brushing'.

Specifications

Model	Sonic style 450	Sonic style 456
Supply voltage	AC220-240V 50Hz	
Power consumption	Approx. 2.5W	
Charging Time/Method	16 hours / Electromagnetic induction	
Time of continuous usage	Under standard swing mode: approx. 20 minutes Automatically stops after 3 minutes	
Weight	Main unit: 62g (include brush head), Charger: 210g	
Modes	Standard Swing- 25,500 strokes/min	Point Swing- 33,000 strokes/min Wide Swing- 13,000 strokes/min Standard Swing- 25,500 strokes/min
Package Contents	Toothbrush main unit, Triple Cleaning Head, charger, instruction manual	Toothbrush main unit, Triple Cleaning Head, Super-Fine Soft Bristle Head, charger, instruction manual
Optional	<ul style="list-style-type: none"> • Duplex Head SB-050 • Triple Cleaning Head SB-070 • Super-Fine Soft Bristle Head SB080 • Point Brush SB-090 	

OMRON

OMRON HEALTHCARE CO., LTD.
24, Yamanouchi Yamanoshita-cho,
Ukyo-ku, Kyoto, 615-0084, Japan

ASIA PACIFIC HQ
OMRON HEALTHCARE SINGAPORE PTE LTD.
83 Clemenceau Avenue #15-05
UE Square Singapore 239920

www.omron-healthcare.com.sg



Sonic style™

SONIC ELECTRIC TOOTHBRUSH

Lightest in Weight

Omron's Unique
Bristle Tip Swing Technology

Effective Removal
of Dental Plaque

A Good Sense of Health

OMRON

ADVANCED TECHNOLOGY FOR OPTIMUM ORAL HYGIENE



Dentists recommend brushing softly by use of brush tip to remove plaque effectively.

Sonic style, powered with Omron's unique Bristle Tip Swing Technology, improves bristle tip movement up to a maximum of 200%*1.

By enhancing the vertical and horizontal bristle tip movement, dental plaque is effectively removed while being gentle on the gums.

*1 Comparison with Standard Swing mode of this product when using the Triple Cleaning Brush (100g load)
*2 Comparison with Duplex Brush

Sonic style™



features designed particularly for daily use


- + Lightweight, only 62g
- + Soft vibration on the hand
- + Ergonomic design for ease of handling
- + Waterproof toothbrush body and charger cover for easy cleaning
- + 3 minute timer automatically switch off after 3 minutes of usage

BRUSH HEAD IMPROVES PLAQUE REMOVAL RATIO BY UP TO 130%*2



Wide Swing – 13,000 strokes/minute

Enhances horizontal movement of bristle tip up to a maximum of 200%*1 to remove dental plaque around periodontal pockets while being gentle on the gums.



Point Swing – 33,000 strokes/minute

Enhances vertical movement of bristle tip up to a maximum of 150%*1 to remove dental plaque on inner side of teeth and gaps between teeth. Removes dental plaque from molars and in-between teeth.



Standard Swing – 25,500 strokes/minute

High-speed strokes to effectively remove plaque.



TRIPLE CLEANING HEAD SB-070

Combination of 3 types of bristles (Spiral, Diamond-Cut bristles and Dual-Merit bristles) to scrape away dental plaque more effectively.



Spiral Bristle

Spirally-arranged bristles can remove plaque even from the side surface of teeth.



Diamond-Cut Bristle

Removes plaque from any angle with 4 edges.




Dual-Merit Bristle

Scrapes plaque out of periodontal pockets.



SUPER-FINE SOFT BRISTLE HEAD SB-080

For soft and gentle brushing, ideal for sensitive gums. With combination of 2 types of bristles (Ultra-Fine and Dual-Merit bristles).



Ultra-Fine Bristle

Reach uneven surface or in-between teeth, with 150 ultra-fine bristles in 1 tuft.



Dual-Merit Bristle


Scrapes plaque out of periodontal pockets.



POINT BRUSH SB-090

Effective interdental brush, for cleaning difficult to reach area.

PLAQUE REMOVING TEST OF BRUSHES



*Source: The Nippon Dental University, Japan (2007)